# beurer

PM 62

## ENGLISH

### Heart Rate Monitor

Instruction manual ..... 1

## ESPAÑOL

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Manual de instrucciones ..... 55

#### **READ THIS MANUAL COMPLETELY AND CAREFULLY BEFORE USING THIS PRODUCT**

Keep this manual in a safe location for future reference

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## 1. IMPORTANT SAFETY NOTES

### Signs and symbols

The following signs appear in the Safety Section (pages 1, 2, 3, 4, 5, 6) and in this manual on pages 11, 15, 16, 22, 27, 33, 36, 40, 43, 44, 45, 46, 47.

**READ THIS ENTIRE MANUAL INCLUDING THE SAFETY SECTION AND ALL INSTRUCTIONS AND WARNINGS COMPLETELY AND CAREFULLY BEFORE USING THIS PRODUCT. FOLLOW ALL SAFETY INSTRUCTIONS AND WARNINGS TO AVOID HAZARDOUS SITUATIONS AND TO MAKE CORRECT USE OF THIS PRODUCT.**

	<b>WARNING</b>	WARNING indicates a hazardous situation which, if not avoided, could result in death or serious injury.
	<b>CAUTION</b>	CAUTION indicates a hazardous situation which, if not avoided, may result in minor or moderate injury.
	<b>NOTICE</b>	NOTICE addresses practices not related to personal injury, such as product and/or property damage.

### **WARNING:**

- This product is a supplementary not a medical device. It is a training instrument, which was developed for the measurement and representation of the human heart rate and should not replace professional medical devices which are used to measure heart rate. The heart rate readings are for reference only, and no responsibility is accepted for the consequences of any erroneous readings.
- The signals used by this product may interfere with a pacemaker or other implanted electronic devices. Consult the manufacturer of the implant device and your physician before using this product. The frequency used by this device

may cause the above mentioned medical devices to malfunction.

- Consult your physician or healthcare provider to help you determine your maximum heart rate, upper and lower limits, the exercise frequency and duration of training appropriate for you and your specific goals. Limits, training schemes, goals, or the like being not in accordance with your health condition represent a risk for your health.
- People suffering from heart and circulatory illnesses should only use this heart rate monitor after having consulted their doctor.
- Always consult your physician or healthcare provider before beginning or modifying a weight reduction or exercise program.
- Always consult your physician to determine what will be suitable for you.
- The fitness test result is for your reference only. Your fitness level can only be ascertained by a licensed physician.
- This product is not designed to detect arrhythmia or irregular rhythms and will interpret them as noise or interference. Arrhythmia can only be ascertained by a licensed physician.
- With high-risk sports, please note that the use of the heart rate monitor can represent an additional source of injury.
- It is important to be sensitive to your body's responses during training. If you feel unexpected pain or excessive fatigue when training, stop the exercise or continue at a lighter intensity. In case of emergency call 911 immediately.
- Keep the unit out of the reach of children, toddlers and infants. This product contains small parts that may present a choking hazard. **NEVER LEAVE CHILDREN OR THOSE WHO REQUIRE CLOSE SUPERVISION UNATTENDED WITH THIS DEVICE.** Immediately contact your physician if any piece is swallowed. In case of emergency call 911 immediately.
- The EasyFit Software is a motivational tool which helps you to monitor your physical activities. Never use it as a basis

for any decision concerning your health, medication, weight reduction or exercise program, or the like. Always consult your physician what will be suitable for you.

- This product is not a toy. Keep it out of reach of children, toddlers, and infants.
- Packaging materials are a deadly hazard for children and can cause suffocation. Remove all packaging materials immediately and keep them away from children at all times.
- Swallowing batteries and/or battery fluid can be extremely dangerous. Keep the battery and unit out of the reach of children and disabled persons. Should any person swallow a battery and/or battery fluid, please call 911 immediately.
- Batteries should not be charged or reactivated by any other means. The batteries may explode.
- Batteries should not be taken apart, thrown in the fire or short circuited. The batteries may explode.
- Should battery fluid leak and come into contact with your eyes or skin immediately rinse with plenty of clean water. Please call 911 immediately.
- The PC interface cable, the chest strap and the elastic tension strap presents a strangulation hazard. Keep this product away from children and those who require close supervision, e.g. people with mental disorders. **NEVER LEAVE CHILDREN OR THOSE WHO REQUIRE CLOSE SUPERVISION UNATTENDED WITH THIS DEVICE.**

## **CAUTION:**

- Use the unit only on human beings and for its intended purpose as outlined in the operating instructions. Do not use for any other purpose. Any form of improper use can be dangerous. The manufacturer is not liable for damages that are caused by improper or incorrect use.
- The device is only intended for private use. Do not use this product for commercial or professional purposes.

- If a battery has leaked do not touch the battery fluid. Avoid skin contact (e.g. put on protective gloves) and clean the battery compartment with a dry cloth.
- Do not disassemble the unit. You could get injured.
- Do not look at the display while you are running to avoid accidents by stumbling.
- Do not adjust the wristband too tight. You could get injured.
- Do not use mobile telephones, microwaves, and other devices which generate strong electrical or electromagnetic fields near the unit. This may result in incorrect operation of the unit.
- Do not use any parts and accessory other than those explicitly recommended by “Beurer North America” for use with this product. Parts and accessories which are not approved for use with the device may cause damage to your health and to the product.

## **NOTICE:**

- Never attempt to repair, open and/or disassemble the unit or adjust it yourself. This may damage the unit and impair the functions. If you need to have the unit repaired, please contact customer service at 1-800-536-0366. Before submitting any complaint, first check the battery of the chest strap and replace it if necessary.
- Changes or modifications to the unit will nullify the user warranty.
- The non-user-replaceable battery of the heart rate monitor has to be changed by a specialist or by Beurer customer service. Improper replacement will damage the unit and void the user warranty.
- Protect the unit against hard knocks (e.g. dropping the unit on the floor), moisture, dirt, dust, chemicals, extreme hot or cold temperatures, major temperature fluctuations, direct exposure to sunlight and heat sources which are too close (e.g. stoves, heating radiators). This may damage the unit

and impair the functions. The device must be stored and used in the specified ambient conditions. Please see chapter "15 Care and Maintenance" and chapter "17 Technical specifications" for details.

- Do not use aggressive solvents, abrasive cleaning agents, detergents or any other strong chemicals to clean the device and accessory. Clean the heart rate monitor with a dry soft cloth.
- Do not jump into water (high pressure load) and press buttons in the rain or while in water. This might impair the impermeability of the monitor to moisture and therefore damage the unit.
- Do not wash the chest strap and the monitor in a washing machine, dish washer or tumble dryer.
- Do not soak, spin-dry, iron, dry clean or bleach the elastic tension strap. Do not use detergent with bleach or fabric softener. This may damage the elastic tension strap.
- Avoid contact with sunscreen or the like since these can cause damage to the imprint or plastic components.
- Do not excessively twist, bend or stretch the chest strap. This may damage the chest strap, its conductive points and impair the functions.
- Do not use any sharp objects to remove the battery from the chest strap.
- Always use the size and type of battery indicated. The use of other batteries will damage the device.
- Never use rechargeable batteries. This may damage the unit.
- Leaking batteries may damage the device. If you do not intend to use the unit for longer periods, remove the battery from the chest strap and have the monitor's battery removed by a specialist before placing the device in storage.
- Do not drop or insert any object into any opening. This may damage the unit

- **Do not press the buttons with excessive force or with any pointed objects.**
- **Observe the local regulations for material disposal. Dispose of the device, components and optional accessories according to applicable local regulations. Unlawful disposal may cause environmental pollution.**
- **Batteries can contain toxins that are harmful to the environment. Always dispose of batteries in accordance with applicable local regulations. Do not dispose of batteries with normal household waste.**

# CONTENT

<b>1. Important Safety Notes .....</b>	<b>1</b>
<b>⚠ WARNING:</b> .....	1
<b>⚠ CAUTION:</b> .....	3
<b>NOTICE:</b> .....	4
<b>2. Package Content .....</b>	<b>10</b>
<b>3. General Information for Training .....</b>	<b>11</b>
<b>4. Functions of the HR monitor .....</b>	<b>13</b>
<b>5. Transmission of signal and methods of Devices Measurement .....</b>	<b>14</b>
5.1 Receiving range of the heart rate monitor .....	14
5.2 Sensors on the chest strap .....	14
5.3 Operating Range .....	14
<b>6. Getting started .....</b>	<b>14</b>
6.1 Activating the HR watch from transport mode .....	14
6.2 Activating heart rate monitor from power saving mode .....	15
6.3 Putting on the heart rate monitor .....	15
6.4 Putting on the chest strap .....	15
6.5 Water Resistance .....	16
<b>7. General operation of the HR monitor .....</b>	<b>16</b>
7.1 Buttons on the HR watch .....	16
7.2 Display .....	17
7.3 Menus .....	18
7.4 Power saving mode .....	18
7.5 Main display .....	18
7.6 Overview of the menus .....	19
<b>8. Basic settings .....</b>	<b>19</b>
8.1 Overview .....	19
8.2 Set units .....	20
8.3 Enter personal data .....	20
8.4 Set training zone .....	22
<b>9. Watch settings .....</b>	<b>24</b>
9.1 Overview .....	24
9.2 Set time of day and date .....	24
9.3 Set alarm clock .....	26
<b>10. Recording of training .....</b>	<b>27</b>
10.1 Overview .....	27
10.2 Stop run times .....	29

10.3 Record laps .....	29
10.4 Stop and store recording .....	30
10.5 Delete all recordings .....	31
10.6 Analyzing memory messages .....	32
<b>11.Results.....</b>	<b>33</b>
11.1 Overview.....	33
11.2 Display training data .....	33
11.3 Lap results .....	35
<b>12.Fitness test.....</b>	<b>36</b>
12.1 Overview.....	36
12.2 Conduct fitness test .....	37
12.3 Display fitness values .....	39
<b>13.EasyFit software .....</b>	<b>40</b>
13.1 Overview.....	40
13.2 Installation of the software .....	41
13.3 Data transfer .....	41
13.4 System requirements .....	42
13.5 EasyFit Instructions .....	42
<b>14.Batteries .....</b>	<b>43</b>
14.1. Heart Rate Monitor .....	44
14.2. Chest Strap.....	44
<b>15.Care and Maintenance .....</b>	<b>45</b>
15.1 Cleaning and Care.....	46
15.2 Storage.....	47
15.3 Disposal .....	47
<b>16.Distortion and interference.....</b>	<b>47</b>
<b>17.Technical specifications .....</b>	<b>48</b>
<b>18.List of abbreviations .....</b>	<b>50</b>
<b>19.Troubleshooting .....</b>	<b>52</b>
<b>20.Warranty .....</b>	<b>53</b>

## **SAVE THESE INSTRUCTIONS FOR FUTURE REFERENCE**

If the instruction manual is damaged or if you no longer have the instruction manual in your possession, please contact Beurer customer service. Please see warranty for service contact.

### **Dear Customer,**

Thank you for choosing this Beurer product. Our name stands for high-quality, thoroughly tested products for applications in the areas of heat, weight, blood pressure, body temperature, pulse, gentle therapy, massage and air.

This heart rate monitor is great tool which can help you to increase the benefit from your exercise program. It helps you to improve the effectiveness of your exercise program by monitoring and quantifying your results every step of the way.

Please read these instructions completely and carefully before using this product and keep them for later reference. Be sure to make them accessible to other users and observe the information they contain.

With best regards,  
Your Beurer Team

## 2. PACKAGE CONTENT

Heart rate monitor with wrist band	
Chest strap	
Elastic tension strap (adjustable)	
Bike mount	
Storage box	
Easy Fit Software and PC interface (USB)	

In addition:

- Operating instructions
- Quick reference guide

**NOTE:** Prior to use, you have to ensure that the heart rate monitor is free of damage. If in doubt, do not use the heart rate monitor and contact Beurer customer service. Please see warranty for service contact.

### 3. GENERAL INFORMATION FOR TRAINING

#### **WARNING:**

- This product is a supplementary not a medical device. It is a training instrument, which was developed for the measurement and representation of the human heart rate and should not replace professional medical devices which are used to measure heart rate. The heart rate readings are for reference only, and no responsibility is accepted for the consequences of any erroneous readings.
- The signals used by this product may interfere with a pacemaker or other implanted electronic devices. Consult the manufacturer of the implant device and your physician before using this product. The frequency used by this device may cause the above mentioned medical devices to malfunction.
- Consult your physician or healthcare provider to help you determine your maximum heart rate, upper and lower limits, the exercise frequency and duration of training appropriate for you and your specific goals. Limits, training schemes, goals, or the like being not in accordance with your health condition represent a risk for your health.
- People suffering from heart and circulatory illnesses should only use this heart rate monitor after having consulted their doctor.
- Always consult your physician or healthcare provider before beginning or modifying a weight reduction or exercise program.
- Always consult your physician to determine what will be suitable for you.
- This product is not designed to detect arrhythmia or irregular rhythms and will interpret them as noise or interference. Arrhythmia can only be ascertained by a licensed physician.
- With high-risk sports, please note that the use of the heart rate monitor can represent an additional source of injury.

This heart rate monitor serves to measure the human pulse. By means of various settings, you can support your individual training program and monitor your pulse. You should monitor your pulse during exercise because you should prevent your heart from excessive strain on the one hand and on the other in order to achieve the optimal training results.

The table provides you with some information on the selection of your training zone. The maximum heart rate reduces with advancing age. Information on the training heart rate must therefore always be in relation to the maximum heart rate. The following rule of thumb helps to identify the maximum heart rate:

### **220 - age = maximum heart rate**

The following example applies to a 40-year-old person:  $220 - 40 = 180$

This rule of thumb offers only an indication, for exact determination of the maximum heart rate consult your doctor or trainer.

	<b>Hearth health programme</b>	<b>Fat burning zone</b>	<b>Fitness programme</b>	<b>Endurance training programme</b>	<b>Anaerobic training</b>
<b>Max. heart rate range</b>	50 – 60%	60 – 70%	70 – 80%	80 – 90%	90 – 100%
<b>Goal</b>	Light exercise	Weight Management	Aerobic base building	Optimal conditioning	Elite athletic training
<b>Possible Effects</b>	Strengthens the heart and the circulatory system	The body burns the highest percentage of calories from fat. Strengthens the heart and the circulatory system, enhances fitness	Improves the respiratory and circulatory system. Ideal to promote the basic endurance	Improves speed maintenance and increases basic speed	Overloads specific muscles. High injury risk for athletes at amateur level, risk related to the heart in case of disease
<b>Suitable for</b>	Ideal for beginners	Weight control and loss	Athletes at amateur level	Ambitious athletes at amateur level, professional athletes	Only professional athletes
<b>Training</b>	Regenerative training		Fitness training	Endurance training	Fitness training addressing specific zones

See the chapter "8 Basic Settings", section "8.4 Set training zone" on how to individually set your training zone. See the chapter "12 Fitness Test" on how the heart rate monitor sets your individual limits of the suggested training zone through the fitness test.

## **4. FUNCTIONS OF THE HR MONITOR**

For product information, questions or customer service call toll-free 1-800-536-0366 or email [info@beurer.com](mailto:info@beurer.com).

### **HR functions**

- Heart rate measurement
- Transmission: analogue
- Individual training zone adjustable
- Acoustic and visual alarm as well as colour change of LED for exit from the training zone
- Average heart rate (AVG)
- Maximum training heart rate
- Calorie consumption in Kcal (basis jogging)
- Fat burning in [oz] / [g]

### **Settings**

- Gender
- Weight
- Unit of weight [kg] / [lb]
- Height
- Unit of height [cm] / [inch]
- Age
- Maximum heart rate

### **Time functions**

- Time of day [12h] / [24h]
- Calendar (up to 2099) and weekday
- Date
- Alarm / alarm clock
- Stop watch
- 50 laps for each training, recall of all of them or each one. Display of lap time and average heart rate

### **Fitness test**

- Fitness index with interpretation
- Maximum oxygen-breath volume (VO<sub>2</sub>max)
- Basal metabolic rate (BMR)
- Active metabolic rate (AMR)
- Calculated maximum heart rate
- Suggested training zones

## 5. TRANSMISSION OF SIGNAL AND METHODS OF DEVICES MEASUREMENT


### 5.1 Receiving range of the heart rate monitor

Your heart rate monitor picks up the heart rate signals sent by the chest strap transmitter within a range of 28 inches (70 cm).

### 5.2 Sensors on the chest strap

The entire chest strap consists of two components: the actual chest strap and an elastic tension strap. On the inside of the chest strap contacting the body, there are two rectangular, corrugated sensors. Both sensors determine your heart rate and transmit this to the heart rate monitor.

**NOTE:** The chest strap is compatible with other analogue devices, using the same transmission frequency, such as e.g. ergometers. Please see chapter "17 Technical specifications" for details.

contact points / sensors

### 5.3 Operating Range

The chest strap and the heart rate monitor should be worn or mounted within 28 inches (70 centimetres) from each other. In case of exceeding this length the monitor displays "0" HR.

## 6. GETTING STARTED

### 6.1 Activating the HR watch from transport mode

During transport the monitor is blank. You can remove the transport mode by pressing the button **MENU** and then **START/STOP**.

## 6.2 Activating heart rate monitor from power saving mode

You can activate your heart rate monitor by pressing any button. If the heart rate monitor has not received any signal for approx. five minutes, it switches to the power saving mode. At that point, only the current time and the date will be displayed. Press any button in order to reactivate the heart rate monitor.

## 6.3 Putting on the heart rate monitor


### CAUTION:

- Do not adjust the wristband too tight. You could get injured.

You can wear your heart rate monitor like a wrist watch. In order to fix the heart rate monitor on bike handle bars, attach the bike mount included to the handle bars. Tighten the heart rate monitor around the mount such that it is not able to slip while riding.

## 6.4 Putting on the chest strap

- Tighten the chest strap with the elastic tension strap. Adjust the length of the elastic tension strap such that the strap is snug, but neither too loose nor too tight. Position the strap around the chest such that the logo is facing outwards and in the correct position directly over the breast bone. **For men**, the chest strap should be located directly below the pectoral muscles, **for women**, directly below the breast.

- Since optimal contact between the skin and the heart rate sensors can not be created immediately, it may take some time until the heart rate is measured and displayed. If necessary, alter the position of the transmitter in order to achieve optimal contact. Most of the time, sufficient contact is achieved through the accumulation of sweat beneath the chest strap. You can, however, also wet the

contact points on the inside of the chest strap. For this, lift the chest strap up from the skin slightly and wet the two sensors with a little saliva or water or a conductive gel.

The contact between the skin and chest strap may not be interrupted also with greater movements of the thorax such as during deep breathing. Significant chest hair can interfere with and even prevent contact.

- Put on the chest strap a few minutes before starting in order to warm it up to body temperature and establish optimal contact.

## 6.5 Water Resistance

### NOTICE:

- **Do not jump into water (high pressure load) and press buttons in the rain or while in water. This might impair the impermeability of the monitor to moisture and therefore damage the unit.**

The PM62 is suitable for bathing and swimming.

Water Resistance = 98 ft. (30 m).

**NOTE:** The consistency of sea water, or chemicals used in a swimming pool may impact the heart rate signal reception.

## 7. GENERAL OPERATION OF THE HR MONITOR

### 7.1 Buttons on the HR watch

1	START/STOP	Starts and stops functions and adjusts values. Press and holding accelerates the entry.
2	OPTION/SET	Switches between the functions within one menu and confirms settings.
3	MENU	Switches between the different menus.
4	Wi-Fi / *	Pressing quickly illuminates the display for five seconds. Pressing and holding activates/ deactivates the acoustic alarm when exiting the training zone.

**NOTE:** If no button is pressed the heart rate monitor will automatically exit the setting mode after 60 seconds and turns back to the main display.

## 7.2 Display

1 Top line

2 Middle line

3 Bottom line

### 4 Symbols

Alarm clock is activated.

Heart rate is being received.

Training zone has been exceeded.

Training zone has been underrun.

Acoustic alarm for exit from the training zone is activated.


The sound when pressing the buttons is activated.

- 5 Two LEDs: Blink in time with heart rate. Colour change from green to red for exit from the training zone.

An explanation of the display abbreviations can be found in the list of abbreviations at the end of the operating instructions.


## 7.3 Menu

You can switch menus with the  Button.

## 7.4 Power saving mode

If your heart rate monitor has not received any signal for approx. five minutes, it switches to the power saving mode. In that case, only the current time and the date will be shown. Press any button to reactivate the HR monitor.


24 h format

12 h format

## 7.5 Main display

When changing menus, the name of the menu is first shown, then your HR monitor will automatically switch to the main display of the active menu. From there, you have to first choose a submenu or you can directly access functions.

## 7.6 Overview of the menus

Menu	Meaning	Description: see chapter...
Time	Time	9 Watch Settings
Training	Training	10 Recording of Training
Result	Result	11 Results
Lap Result	Lap result	11.3 Lap results
Setting	Settings	8 Basic settings
FitTest	Fitness test	12 Fitness test

## 8. BASIC SETTINGS

### 8.1 Overview

In the menu **Setting** you can:

- Change the measurement units for height (cm/inch) and weight (kg/lb).
- Enter your personal data (user). Based on this data, your HR monitor determines your calorie consumption and your fat burning during training!
- Set your training zone (limits). See chapter "3 General Information for Training" for explanations on how to set the training zones. Your HR monitor issues an acoustic and visual alarm during your training when the training zone is exited.

Once you are in the menu **Setting**, your HR monitor automatically switches to the main display. This displays **User**.


138

Setting

138



User

## 8.2 Set units

Switch to the menu <b>Setting</b> by pressing the <b>MENU</b> button. Wait till the monitor automatically switches to the main display. This displays <b>User</b> .	
Press <b>OPTION/SET</b> button. <b>Limits</b> is displayed.	
Press <b>OPTION/SET</b> button. <b>Units</b> is displayed.	
Press <b>START/STOP</b> button. <b>kg/lb</b> is displayed. The weight unit blinks. Press <b>START/STOP</b> button to switch between the units. Confirm with <b>OPTION/SET</b> button.	
<b>cm/inch</b> is displayed. The unit of measurement for height blinks. Press <b>START/STOP</b> button to switch between the units. Confirm with <b>OPTION/SET</b> button. <b>Units</b> is displayed.	

## 8.3 Enter personal data

Switch to the menu <b>Setting</b> by pressing the <b>MENU</b> button. Wait till the monitor automatically switches to the main display. This displays <b>User</b> .	
Press <b>START/STOP</b> button. <b>Gender</b> is displayed. Press <b>START/STOP</b> button to set your gender. Select <b>M</b> , if you are male. Select <b>F</b> , if you are female. Confirm with <b>OPTION/SET</b> button.	

<p><b>AGE</b> is displayed.</p> <p>Set your age (10-99) by pressing the <b>START/STOP</b> button. Confirm with <b>OPTION/SET</b> button.</p>	
<p><b>MaxHR</b> is displayed. The HR monitor automatically calculates the maximum heart rate according to the formula: <b>220 - age = MaxHR</b></p> <p><b>Please note</b> that the HR monitor also determines the MaxHR after every fitness test and sets it as a suggested value.</p> <p>You can overwrite the value calculated by pressing the <b>START/STOP</b> button (setting range: 65 to 239). Confirm with <b>OPTION/SET</b> button.</p>	
<p><b>Weight</b> is displayed.</p> <p>Press <b>START/STOP</b> button to select your weight. Setting range: 44 to 499 lbs (20 to 226 kg). Confirm with <b>OPTION/SET</b> button.</p>	
<p><b>Height</b> is displayed.</p> <p>Press <b>START/STOP</b> button to set your height. Setting range: 2'6" to 7'6" feet (80 to 227 cm). Confirm with <b>OPTION/SET</b> button. <b>User</b> is displayed.</p>	

## 8.4 Set training zone

### **WARNING:**

- **Consult your physician or healthcare provider to help you determine your maximum heart rate, upper and lower limits, the exercise frequency and duration of training appropriate for you and your specific goals. Limits, training schemes, goals, or the like being not in accordance with your health condition represent a risk for your health.**

You can set your individual upper and lower heart rate limit as described below.

**NOTE:** Both limits are automatically set as a training zone suggestion after every fitness test and that the old data are consequently overwritten. For this, the lower limit of the suggested training zone totals 65%, the upper limit 85% of the maximum heart rate.

You can calculate the lower and upper heart rate limits by multiplying your maximum heart rate by the percentages for the selected training zone (see chapter "3 General Information for Training"). Please note that the lower limit must be smaller than the upper limit.



E.g. a 50-year-old man training for Weight Management (fat burning zone).

His lower heart rate limit:  $(220-50)*65\% = 110.5 \rightarrow 111$

His upper heart rate limit:  $(220-50)*85\% = 144.5 \rightarrow 145$

**NOTE:** For a safe and efficient training determining the individual training zone is a critical step.

Do not precipitate – Start slowly. Gradually step up your training zone as you become more physically fit. Select the training zone which suits your needs the most.

<p>Switch to the menu <b>Setting</b> by pressing the <b>MENU</b> button. Wait till the monitor automatically switches to the main display. This displays <b>User</b>.</p>	
<p>Press <b>OPTION/SET</b> button. <b>Limits</b> is displayed.</p>	
<p>Press <b>START/STOP</b> button. <b>On LED</b> is displayed. The LED is activated (ON) or shut off (OFF) with pressing <b>START/STOP</b> button.</p> <ul style="list-style-type: none"> <li>■ <b>LED (ON)</b> The LED gives a red light when your heart rate is outside the training zone. The LED gives a green light when your heart rate is inside the training zone.</li> <li>■ <b>LED (OFF)</b> LED is deactivated. No light will flash up.</li> </ul> <p>Confirm with <b>OPTION/SET</b> button.</p>	
<p><b>Lo limit</b> is displayed. Set the lower limit for your training zone by pressing <b>START/STOP</b> button. Setting range: 40-238 Confirm with <b>OPTION/SET</b> button.</p>	
<p><b>Hi limit</b> is displayed. Set the upper limit for your training zone by pressing <b>START/STOP</b> button. Setting range: 41-239 Confirm with <b>OPTION/SET</b> button and to return to the display <b>Limits</b>.</p>	

## 9. WATCH SETTINGS

### 9.1 Overview

In the **Time** menu you can:

- Set the time and the date.
- Set the alarm clock.

When you are in the **Time** menu the heart rate monitor automatically activates the main display. In this mode both the weekday and the current date are displayed.

### 9.2 Set time of day and date

Switch to the menu **Time** by pressing the **MENU** button.


#### Set time format

Press **START/STOP** button for approx. 5 seconds. **24 hr** or **12 hr TimeFORM** is displayed. Press **START/STOP** button to choose the 12-hour AM/PM or 24-hour universal clock. Confirm with **OPTION/SET** button.

#### NOTE:


- If you select 24hr time format date will be displayed “weekday day.month”.
- If you select 12 hr time format date will be displayed “weekday month.day”.

<p><b>Hour</b> is displayed. The hour setting blinks. Set the hour with the <b>START/STOP</b> button.</p> <p>Confirm with <b>OPTION/SET</b> button.</p>	
<p><b>Minute</b> is displayed. The minute setting blinks. Press <b>START/STOP</b> button to set the minute (0-59).</p> <p>Confirm with <b>OPTION/SET</b> button.</p>	
<p><b>Year</b> is displayed. The year setting blinks. Press <b>START/STOP</b> button to set the year. The calendar goes up to 2099.</p> <p>Confirm with <b>OPTION/SET</b> button.</p>	
<p><b>Month</b> is displayed. The monthly setting blinks. Press <b>START/STOP</b> button to set the month (1-12).</p> <p>Confirm with <b>OPTION/SET</b> button.</p>	
<p><b>Day</b> is displayed. The day setting blinks. Press <b>START/STOP</b> button to set the day (1-31).</p> <p>Confirm with <b>OPTION/SET</b> button.</p> <p>The display shows the current time and date.</p>	

### 9.3 Set alarm clock

When the alarm clock function is activated, the alarm will go off at the set time every day. You can turn off the alarm clock by pressing any button.

<p>Switch to the menu <b>Time</b> by pressing the <b>MENU</b> button.</p>	
<p>Press <b>OPTION/SET</b> button. <b>Alarm</b> is displayed.</p>	
<p>Press <b>START/STOP</b> button for approx. 5 seconds. <b>On Alarm</b> is displayed. Press <b>START/STOP</b> button to turn the alarm clock off (Off) and on (On). Confirm with <b>OPTION/SET</b> button.</p>	
<p><b>Hour</b> is displayed. The hour setting blinks. Press <b>START/STOP</b> button to set the hour. Confirm with <b>OPTION/SET</b> button.</p>	
<p><b>Minute</b> is displayed. The hour setting blinks. Press <b>START/STOP</b> button to set the minute. Confirm with <b>OPTION/SET</b> button and to return back to the main display.</p>	

## 10. RECORDING OF TRAINING

### **WARNING:**

- The signals used by this product may interfere with a pace-maker or other implanted electronic devices. Consult the manufacturer of the implant device and your physician before using this product. The frequency used by this device may cause the above mentioned medical devices to malfunction.
- It is important to be sensitive to your body's responses during training. If you feel unexpected pain or excessive fatigue when training, stop the exercise or continue at a lighter intensity. If case of emergency call 911 immediately.

### **CAUTION:**

- Do not look at the display while you are running to avoid accidents by stumbling.
- Do not use any parts and accessory other than those explicitly recommended by "Beurer North America" for use with this product. Parts and accessory which are not approved for use with the device may cause damage to your health and to the product.

### 10.1 Overview

#### **NOTE:**

- Before starting training recording, you need to wear the chest strap.
- Avoid doing any awkward motion during measurement. It will create inaccurate heart rate reading.
- In addition to training intensity, medication for heart conditions / blood pressure / physiological conditions / asthma / breathing, etc. as well as some energy drinks, alcohol, and nicotine may also affect heart rate.

In the menu **Training** you can:

- Stop your run times and thereby record your training session.
- Record individual laps.
- Stop and save your recording.
- Delete all recordings.
- Analyze memory messages.

Once you are in the menu **Training**, your HR monitor automatically switches to the main display. This displays the already recorded training time on the middle line (after Save/Delete: 00:00:00) and on the line below **Start**.

## Things to know about the recordings

- As soon as you have started the recording, your HR monitor starts recording your training data. Your HR monitor continually records your training data until you save or delete this. This also applies when you first start training on the following day. The running time also allows you to stop your training session at any time.
- Your HR watch has a memory. The recording time is up to 24 hours and depends on the number of laps. A maximum of 50 laps is possible for each training. Your HR watch informs you if the memory is almost full. At this point, at the latest, we recommend transmitting the recordings to your PC and thus emptying the memory on the HR watch. Alternatively you can delete all recordings manually. Read the section "10.5 Delete all recordings" in this chapter to find out how to delete the recordings.

- You can let your training data be displayed. The last recording saved, can be recalled in the menu **Result**. In the menu **Lap Result** you can see the laps for this recording.
- As soon as you have transferred the recordings into the EasyFit software, you can view and analyze them there. The memory will be deleted every time you have successfully transferred the data to the PC. You can find explanations regarding the transfer of data to the PC in the chapter "13 EasyFit Software".

## 10.2 Stop run times

<p>Switch to the menu <b>Training</b> by pressing the  button. If you have already stopped a run time, the recorded training time is displayed on the middle line.</p>	
<p>Start the recording by pressing the  button.</p> <p><b>Run</b> is displayed. Even if you change to another menu during the recording, the time continues to run in the background.</p>	
<p>Stop the recording by pressing the  button. <b>Stop</b> is displayed. You can interrupt and resume your recordings at any time by pressing the  button.</p>	

## 10.3 Record laps


While the recording of your entire training is in progress, you can record individual laps. The time continues to run as the laps are recorded.

<p>You are in the menu <b>Training</b> and your training is currently being recorded. <b>Run</b> is displayed.</p>	
<p>Press <b>OPTION/SET</b> button to time a new lap. The display shows the time and number of the last lap timed for five seconds, for example <b>Lap 1</b>.</p>	
<p>The display automatically switches back to <b>Run</b> and shows the number of the current lap, for example <b>Run 2</b>.</p>	
<p>Press <b>OPTION/SET</b> button to record additional training laps as needed. You can record up to a total of 50 laps for each training.</p>	

## 10.4 Stop and store recording

Please note that after finishing and saving, you can no longer have the data displayed in the **Result** and **Lap Result** menus. You can only display them on your PC with the help of the EasyFit Software.


<p>You are in the menu <b>Training</b>. You have just recorded a training. <b>Stop</b> is displayed.</p>	
--	--

<p>Save the recording by pressing <b>START/STOP</b> button for approx. 5 seconds.</p> <p><b>Hold Save</b> is displayed.</p>	
<p>The display automatically switches back to the main display.</p> <p>The data is now ready to be transmitted to the PC.</p>	

## 10.5 Delete all recordings


Your HR watch informs you automatically when the memory is almost full. Before you then delete the memory, we recommend analyzing the training recording or transmit to the PC! The memory on the HR watch will be deleted automatically every time you successfully transfer data to the PC. Alternatively you can delete the memory manually as described here.

The automatic or manual deletion causes all training recordings including your laps to be deleted.

<p>You can delete the data <b>before</b> or <b>after</b> the training:</p> <ul style="list-style-type: none"> <li>■ You are in the menu <b>Training</b>. <b>Start</b> is displayed.</li> <li>■ You are in the menu <b>Training</b>. You have just recorded a training. <b>Stop</b> is displayed.</li> </ul>	
<p>Delete all recordings by pressing <b>OPTION/SET</b> button for approx. 5 seconds.</p> <p><b>MEM Delete</b> is displayed.</p>	

The monitor automatically switches to the main display.

138  
0000.00  
Start

## 10.6 Analyzing memory messages

In the menu **Training** your HR watch automatically informs you if the memory is almost full or if almost all the laps are occupied.

Message	Meaning
	When switching to the menu <b>Training</b> by pressing the <b>MENU</b> button the HR watch signals <b>MEM LOW</b> , if less than 10% of the memory or less than 10 laps are free.
	During training, the HR watch signals every 5 seconds <b>MEM x% Free</b> , if less than 10% of the memory is free. E. g. on the display: <b>MEM 8% Free</b> (still 8% free).
	During training the HR watch signals every 5 seconds <b>LAP Rest x</b> , if there are less than 10 (of a total of 100) laps . E. g. on the display: <b>LAP Rest 7</b> (still 7 laps free).

If the memory is full, you cannot start a new training session. Read in the previous section "10.5 Delete all recordings" to find out how to delete the memory manually.

The memory is also deleted every time the data is transferred to the PC. You can find the relevant explanations in the chapter "13 EasyFit Software".

## 11. RESULTS

### WARNING:

- This product is a supplementary not a medical device. It is a training instrument, which was developed for the measurement and representation of the human heart rate and should not replace professional medical devices which are used to measure heart rate. The heart rate readings are for reference only, and no responsibility is accepted for the consequences of any erroneous readings.
- This product is not designed to detect arrhythmia or irregular rhythms and will interpret them as noise or interference. Arrhythmia can only be ascertained by a licensed physician.

### 11.1 Overview

In the menu **Result** you can display your training data of the last recording stopped: mean and maximum heart rate, times within and outside of the training zone as well as calorie consumption and fat burning.

Result

Once you are in the menu **Result**, your HR monitor automatically switches to the main display.

This displays **bPM AVG**.


AVG

### 11.2 Display training data

Switch to the menu **Result** by pressing the

MENU

button.

<p><b>BPM AVG</b> (average training heart rate) is displayed.</p>	
<p>Press the <b>OPTION/SET</b> button.</p> <p><b>HRmax</b> (maximum training heart rate) is displayed.</p>	
<p>Press the <b>OPTION/SET</b> button.</p> <p><b>In [min]</b> (time within the training zone) is displayed.</p>	
<p>Press the <b>OPTION/SET</b> button.</p> <p><b>Lo [min]</b> (time below the training zone) is displayed.</p>	
<p>Press the <b>OPTION/SET</b> button.</p> <p><b>Hi [min]</b> (time above the training zone) is displayed.</p>	
<p>Press the <b>OPTION/SET</b> button.</p> <p><b>[kcal]</b> (the total calorie consumption during the recording of the training session) is displayed. The calculation is conducted based on moderately fast jogging. The basis for the display of the calorie consumption and the fat burning are:</p>	




<ul style="list-style-type: none"> <li>■ Personal data (basic settings)</li> <li>■ the time and heart rate measured during the training session.</li> </ul>	
<p>Press the <b>OPTION/SET</b> button.</p> <p><b>Fat[g]/[oz]</b> (Fat burning during the entire training recording) is displayed.</p>	
<p>Press the <b>OPTION/SET</b> button.</p> <p>Back to the main display.</p>	

## 11.3 Lap results

### Overview

<p>In the menu <b>LAP Result</b> you can have the time and average heart rate as well as the last stopped training recording and the individual laps displayed.</p>	
<p>Once you are in the menu <b>LAP Result</b>, your HR monitor automatically switches to the main display. The top line shows the average heart rate, the middle line the training time last recorded.</p>	

## Display lap results

Switch to the menu <b>LAP result</b> by pressing the  button.	
Press the  button. <b>AVG 1:</b> Results in lap 1 are displayed: The top line shows the average heart rate, the middle line displays the lap time.	
Pressing  button again switches to the individual result display for the following lap. Maximum 50 laps for each training.	
Press the  button to get back to the main display.	

## 12. FITNESS TEST

### **WARNING:**

- The fitness test result is for your reference only. Your fitness level can only be ascertained by a licensed physician.

### 12.1 Overview

In the menu <b>FitTest</b> you can: <ul style="list-style-type: none"><li>■ Test your personal fitness in a 1 mile (1600 meters) run.</li><li>■ Have your fitness evaluated: Fitness index (1-5), maximum oxygen-breath volume (VO<sub>2</sub>max), basal metabolic rate (BMR) and active metabolic rate (AMR).</li></ul>	
---	--


- Have your individual upper and lower limit calculated automatically as a recommended training zone and overwritten the two limits in the menu **Setting**. Also see the section "8.4 Set training zone".
- Have your maximum heart rate (MaxHR) automatically calculated and have it overwritten in the menu **Setting**.

The fitness test is not just designed for active athletes, but also for those with a lower fitness level from all age groups.

One second after entering the Fitness test mode, the main display appears. The top line displays your fitness index:

- 0 = No test has been conducted yet
- 1 = Poor
- 2 = Fair
- 3 = Average
- 4 = Good
- 5 = Elite


On the middle line, the date of your last fitness test is displayed. If you have not yet completed a fitness test, then the preset date is displayed.

## 12.2 Conduct fitness test

The units displayed [mile] or [m] (meter) result from the basic settings. In the chapter "8 Basic Settings", section "8.2 Set units", you can read how to change these units.

Put on the chest strap and wait till the monitor receives a heart rate. Switch to the menu **Fit-Test** by pressing the **MENU** button.

The scrolling display appears **PressSTART to begin**. Your pulse will now be shown on the display and you have five minutes to warm up.

Your HR monitor signals the end of the 5-minute warmup phase with a beep. In the event you would like to shorten the warm-up phase, you can press the **START/STOP** button to switch to the next screen.

Briskly walk exactly 1 mile (1600m). For this, choose a level route of exactly 1 mile (1600m).

The scrolling display **PressSTOP after 1mile (1600 m)** reminds you to press the **START/STOP** button after 1 mile (1600m).

The scrolling display **Test Finished** appears and the time required is displayed.



Press **START/STOP** button to get back to the main display. It displays your fitness index (1-5) with evaluation and the current date.

## 12.3 Display fitness values

The display of the fitness values only makes sense, if you have completed at least one complete fitness test. The fitness values displayed refer to the last test.

<p>Change to the menu <b>FitTest</b> by pressing the <b>MENU</b> button. The fitness index (1-5) and the date of your last fitness test is shown. If no fitness test has been completed yet, a <b>0</b> is displayed.</p>	
<p>Press the <b>OPTION/SET</b> button. <b>VO<sub>2</sub>max</b> is displayed. The maximum oxygen-breath volume is displayed in millilitres within one minute per 2.2 pounds (1 kg) of body weight. It is the basis for the evaluation of your fitness test.</p>	
<p>Press the <b>OPTION/SET</b> button.</p> <p><b>BMR</b> (Basal Metabolic Rate) is displayed. The basal metabolic rate is calculated based on the values of age, gender, height and weight. The BMR specifies the calories needed to maintain your vital functions - in the absence of movement.</p>	
<p>Press the <b>OPTION/SET</b> button.</p> <p><b>AMR</b> (Active Metabolic Rate) is displayed. The active metabolic rate is calculated based on the basic metabolic rate and an "activity premium", which is dependant on the fitness level.</p>	
<p>Press the <b>OPTION/SET</b> button to get back to the main display.</p>	

## 13. EASYFIT SOFTWARE

### **WARNING:**

- The EasyFit Software is a motivational tool which helps you to monitor your physical activities. Never use it as a basis for any decision concerning your health, medication, weight reduction or exercise program, or the like. Always consult your physician what will be suitable for you.
- The PC interface cable, the chest strap and the elastic tension strap presents a strangulation hazard. Keep this product away from children and those who require close supervision, e.g. people with mental disorders. **NEVER LEAVE CHILDREN OR THOSE WHO REQUIRE CLOSE SUPERVISION UNATTENDED WITH THIS DEVICE.**

### 13.1 Overview

With the EasyFit software you can transfer your training data to a PC via your Heart Rate Monitor's PC interface. The required accessories, consisting of a CD with the PC software EasyFit and a communication cable, are included in the scope of supply. The EasyFit software helps you to monitor your training data and control your achievements over a longer period of time. In addition to that, you will find many other functions for the support of your weight management and workout planning.

- Calendar functions
- Weekly workout planning
- Various analyses and graphics for observing the training progress

**NOTE:** The EasyFit Software can also be used with EasyFit ready Beurer Activity sensor and scales. You thus get an overview of your activity, training and weight data. Please contact Beurer customer service for more information. Please see warranty for service contact.

## 13.2 Installation of the software

To install the software, please pay attention to printed instructions or to the file "Install\_Readme" on the enclosed CD. Please follow the instructions carefully to install the software exactly and in the order specified.

1. Place the CD-ROM supplied in the drive and choose the CD drive in Windows Explorer or on the desktop of your PC.
2. Double-click the "setup.exe" file to run it.
3. Choose your language and follow the installation instructions. Confirm with "Next" or "Install" as appropriate.
4. When you click "Finish", the software starts.
5. Enter your personal settings and details under "New registration" in the EasyFit program and click "Save" or "Next" to confirm.

## 13.3 Data transfer

1. Start EasyFit.
2. Create a new person or click on an existing account.
3. Go to "Data Manager" and click the "Transfer data" button.
4. Connect the PC clip to the heart rate monitor and after that connect the USB plug to the USB port of your PC.
5. After successful data transfer your training data can be viewed in the menu "Training", "Data Manager" or can be called up individually in the calendar.

**NOTE:** Do not turn off the PC or unplug the USB connector during uploading of the measured data.

## FURTHER INFORMATION

Should you have any questions concerning the use of the EasyFit Software please call customer service at toll-free **1-800-536-0366** or email **info@beurer.com**.

## 13.4 System requirements

Operating system	Windows 2000/XP (SP2), Vista, Windows 7
CPU	optimal at 1.2 GHz+
Main memory	optimal at 256 MB RAM+
Hard disk space	optimal 200 MB+ free space – at least 100 MB
Graphic resolution	starting from 1024 x 768
Other requirements	CD-ROM drive, USB interface

## 13.5 EasyFit Instructions

Detailed instructions concerning the different functions and the operation of the software are included in the help function.

## 14. BATTERIES

### **WARNING:**

- Swallowing batteries and/or battery fluid can be extremely dangerous. Keep the battery and unit out of the reach of children and disabled persons. Should any person swallow a battery and/or battery fluid, please call 911 immediately.
- Batteries should not be charged or reactivated by any other means. The batteries may explode.
- Batteries should not be taken apart, thrown in the fire or short circuited. The batteries may explode.
- Should battery fluid leak and come into contact with your eyes or skin immediately rinse with plenty of clean water. Please call 911 immediately.

### **CAUTION:**

- If a battery has leaked do not touch the battery fluid. Avoid skin contact (e.g. put on protective gloves) and clean the battery compartment with a dry cloth.

### **NOTICE:**

- The non-user-replaceable battery of the heart rate monitor has to be changed by a specialist or by Beurer customer service. Improper replacement will damage the unit and void the user warranty.
- Do not use any sharp objects to remove the battery from the chest strap.
- Always use the size and type of battery indicated. The use of other batteries will damage the device.
- Never use rechargeable batteries. This may damage the unit.
- Leaking batteries may damage the device. If you do not intend to use the unit for longer periods, remove the battery from the chest strap and have the monitor's battery removed by a specialist before placing the device in storage.

- **Batteries can contain toxins that are harmful to the environment. Always dispose of batteries in accordance with applicable local regulations. Do not dispose of batteries with normal household waste.**

**NOTE:** This product contains a CR Type Coin Lithium Battery which contains Perchlorate Material – special handling may apply.

See: [www.dtsc.ca.gov/hazardouswaste/perchlorate](http://www.dtsc.ca.gov/hazardouswaste/perchlorate).



### 14.1. Heart Rate Monitor

If the intensity of the display reduces or no display appears at all, have the battery of the HR monitor only changed by a specialist (e.g. clock-maker's shop or jewellery store) or by Beurer customer service otherwise the water resistance might be impaired. Improper replacement will void the warranty.

**NOTE:** Excessive use of the LED light drains the battery more rapidly.

### 14.2. Chest Strap

If the heart rate signals become unsteady or stop your chest strap battery may need to be replaced.

1	You can replace the battery in the chest strap on your own. The battery compartment is located on the inner side. Use a coin to open the compartment.	
2	The positive pole of the battery inserted must face upwards.	
3	Carefully close the battery compartment again.	

## 15. CARE AND MAINTENANCE

### **CAUTION:**

- Do not disassemble the unit. You could get injured.
- Do not use any parts and accessory other than those explicitly recommended by "Beurer North America" for use with this product. Parts and accessory which are not approved for use with the device may cause damage to your health and to the product.

### **NOTICE:**

- Never attempt to repair, open and/or disassemble the unit or adjust it yourself. This may damage the unit and impair the functions. If you need to have the unit repaired, please contact customer service at 1-800-536-0366. Before submitting any complaint, first check the battery and replace them if necessary.
- Changes or modifications to the unit will nullify the user warranty.
- Protect the unit against hard knocks (e.g. dropping the unit on the floor), moisture, dirt, dust, chemicals, extreme hot or cold temperatures, major temperature fluctuations, direct exposure to sunlight and heat sources which are too close (e.g. stoves, heating radiators). This may damage the unit. The device must be stored and used in the specified ambient conditions. Please see chapter "15 Care and Maintenance" and chapter "17 Technical specifications" for details.
- Do not use aggressive solvents, abrasive cleaning agents, detergents or any other strong chemicals to clean the device and accessory. Clean the heart rate monitor with a dry soft cloth.
- Do not jump into water (high pressure load) and press buttons in the rain or while in water. This might impair the im-

permeability of the monitor to moisture and therefore damage the unit.

- Do not wash the chest strap and the monitor in a washing machine, dish washer or tumble dryer.
- Do not soak, spin-dry, iron, dry clean or bleach the elastic tension strap. Do not use detergent with bleach or fabric softener. This may damage the elastic tension strap.
- Avoid contact with sunscreen or the like since these can cause damage to the imprint or plastic components.
- Do not excessively twist, bend or stretch the chest strap. This may damage the chest strap, its conductive points and impair the functions.
- Do not drop or insert any object into any opening. This may damage the unit
- Do not press the buttons with excessive force or with any pointed objects.

## 15.1 Cleaning and Care

From time to time, carefully clean your **chest strap, elastic tension strap and heart rate monitor** with a soap and water solution. Then rinse off all components with clean water. Carefully dry them only with a soft towel. Always thoroughly dry off any moisture on the components before put them into storage.

**NOTE:** Dirt impairs the elasticity and functioning of the transmitter. To avoid build up of material that can interfere with the chest strap's ability to receive a valid pulse, rinse the chest strap under running water after every use. Carefully dry it only with a soft towel.

You can wash the **elastic tension strap** in the washing machine at 30° without fabric softener. Use a washing pouch. The elastic tension strap is not made for the dryer! Detach the chest strap from the elastic tension strap. Never put the **chest strap** in the washing machine or dryer.

## 15.2 Storage

Store your monitor and accessories in the original storage box in a cool and dry place.

**NOTE:** Never store the **chest strap and elastic tension strap** when they are wet and make sure that the contact points of the chest strap do not contact anything damp, such as sport towel or wet elastic tension strap. Sweat and moisture can keep the contact points wet and consequently activate the chest strap, which shortens battery life.

## 15.3 Disposal

### **NOTICE:**

- Observe the local regulations for material disposal. Dispose of the device, components and optional accessories according to applicable local regulations. Unlawful disposal may cause environmental pollution.
- Batteries can contain toxins that are harmful to the environment. Always dispose of batteries in accordance with applicable local regulations. Do not dispose of batteries with normal household waste.

## 16. DISTORTION AND INTERFERENCE

### **WARNING:**

- The signals used by this product may interfere with a pacemaker or other implanted electronic devices. Consult the manufacturer of the implant device and your physician before using this product. The frequency used by this device may cause the above mentioned medical devices to malfunction.

In the HR monitor, a receiver for signals from the chest strap transmitter is integrated. During use, it may arise that the HR monitor receives signals other than those transmitted by the chest strap.

## How can I recognize interference?

Interference can be recognized in that the values displayed are unrealistic. If your normal heart rate amounts to e.g. 130 bpm and it is suddenly shown to be 200, then a short-term disturbance is present. Once you have removed yourself from the vicinity of the source causing the disturbance, the values will automatically resume being accurate. Please note that distorted values can influence the training result displays and overall values.

## What kind of causes for interference are there?

In some situations, short-range electromagnetic interference can occur. In order to avoid these, keep well away from these sources.

- Nearby power lines
- Walking through security gates in airports.
- Nearby overhead lines (tram)
- Nearby traffic lights, television sets, car motors, bike computers
- Nearby exercise equipment
- Nearby mobile phones

## Interference from other devices for measuring heart rates

The HR monitor as a receiver can react to interfering stray signals such as other heart rate monitors within a vicinity of one metre. This can result in incorrect measurements particularly when training in groups. If you e.g. jog next to one another, it can be beneficial to wear the HR monitor on the other wrist or to increase the distance to one's training partner.

## 17. TECHNICAL SPECIFICATIONS

Model	PM62
-------	------

### Heart rate monitor

Real time clock	12/24 hour with hour/minute/second
Heart rate measurement	35 to 240 bpm
Accuracy of heart rate	+/- 1bpm, whichever larger, definition applies to steady state conditions

Water Resistant (not actually keys)	98 ft. (30 m) suitable for swimming and bathing <b>NOTE:</b> The impermeability of the monitor to moisture cannot be guaranteed if you jump into water (high pressure load) and also if the buttons are pressed in the rain.
Battery HR monitor (non-user-replaceable battery)	3V lithium battery, type CR2032 (Recommendation: SONY batteries). <b>NOTE:</b> The service life is approx 15 months if the monitor with pulse display is used for 1 hour daily. In addition, activated functions such as e.g. acoustic alarm or activated LED display could reduce the service life of the battery to up to seven months. The included batteries are only for test mode and could be consumed within one year.

### Chest strap

Chest strap size	max. 51.2" (max. 130 cm)
Water Resistant	Waterproof (suitable for swimming)
Transmission frequency	5.3 KHz
Battery chest strap	3V lithium battery, type CR2032 Service life: approx 28 months

Weight	HR monitor (including battery) = 1.7 oz (47 g) Chest strap (including battery) = 1.6 oz (46 g) Elastic tension strap = 0.9 oz (26 g)
Operating Temperature	14°F to 140° F (-10° C to 60 °C)
Storage and transport temperature	-4°F to 158° F (-20° C to 70 °C)

Subject to technical modifications.

## 18. LIST OF ABBREVIATIONS

Abbreviation	English
12hr	12-hour AM/PM clock
24hr	24-hour universal clock
AGE	Age
Alarm	Alarm
AM	Ante meridiem: time between midnight and noon
AMR	Calories needed for the activity
AVG	Average heart rate during training
AVG Run	Average heart rate during training
BMR	Basal Metabolism
BPM	Heart beats per minute
bPM AVG	Average heart rate during training, beats per minute
bPM Hrmax	Maximal heart rate during training, beats per minute
cm/inch	Centimeter/Inch
Day	Day
f	female
Fat [oz]/[g]	Fat in ounces / gram
FitTest	Fitness Test
Fri	Friday
Gender	Gender
Height	Height
Hi	Time above the training zone
Hi limit	Highest training limit
Hold Save	Save the recording
Hour	Hour
HR	Heart Rate
HRmax	Max. heart rate during training
in	inch
In	Time within the training zone
[Kcal]	Calorie consumption in Kcal
kg	Kilogram
kg/lb	Kilogram/pound

<b>Abbreviation</b>	<b>English</b>
lb	Pound
lap	Lap
LED	Light emitting diode
Limits	Training limits
Lo	Time below the training zone
Lo limit	Lowest training limit
m	male
MaxHR	Maximum Heart Rate (formula: 220 – age)
MEM	Memory
MEM Delete	Memory Delete
Minute	Minute
[min]	Minute
Mon	Monday
Month	Month
PM	Post meridiem: time between noon and midnight
Result	Result
Run	Run time is measured
Sat	Saturday
Setting	Settings (Menu)
SEt	Setting
Start	Start
Stop	Stop
Sun	Sunday
Time	Time (Menu)
Thu	Thursday
Training	Training (Menu)
Tue	Tuesday
Units	Units
User	User
VO <sub>2</sub> max	Maximal Amount of Oxygen
Wed	Wednesday
Weight	Weight
Year	Year

## 19. TROUBLESHOOTING

Problem/Question	Possible Cause / Remedy
No heart rate signal	Shorten the distance between the watch and the chest strap.
	Adjust the position of the chest strap and the elastic tension strap. Make sure the chest strap sensors are wet.
	Check the battery in the chest strap and replace it if necessary.
	Move away from other heart rate monitors or other sources of electrical or radio interference.
	Tighten the elastic tension strap. If the chest strap is too loose, the movement of the contact points disturbs the detection of the heart rate signal.
	Hairy chest may also weaken the contact. Try to find the best possible position for the contact points.
	Dirty contact points. Clean the chest strap. See section "15.1 Cleaning and Care".
Training cannot be started.	Full memory. Delete your saved training sessions or transmit them to the PC.
Unusual high or low values	Move away from other heart rate monitors or other sources of electrical or radio interference.
	Wear the unit properly according to the manual.
	Battery of the chest strap is low/empty.
Obscure display	Low battery. Have the battery changed by a specialist or Beurer customer service.
	Ambient temperature is too low. The display will recover to its normal status when the temperature rises.
LCD is black	The unit has been exposed too long to direct sunlight or high temperatures. It will recover to its normal status under the shade.

<b>Problem/Question</b>	<b>Possible Cause / Remedy</b>
No display	Battery is empty. Have the battery changed by a specialist or Beurer customer service. The unit is broken. Contact customer service at 1-800-536-0366 or email <a href="mailto:info@beurer.com">info@beurer.com</a> .
Mem 0%	Memory is full. As a result the training cannot be started / reset to 00:00. Transfer data to the PC (see chapter "13 EasyFit Software") or delete the memory (see section "10.5 Delete all recordings").

## **20. WARRANTY**

### **Limited Lifetime Warranty For Original Purchaser**

Your Beurer Heart Rate Monitor, Model PM62, excluding the batteries and all accessories, is warranted to be free from defects in materials and workmanship for the life of the product under normal conditions of intended use and service. This warranty extends only to the original retail purchaser and does not extend to retailers or subsequent owners.

We will, at our option, repair or replace the Beurer Heart Rate Monitor, Model PM62, without additional charge, for any part or parts covered by these written warranties. No refunds will be given. Repair or replacement is our only responsibility and your only remedy under this written warranty. If replacement parts for defective materials are not available, Beurer reserves the right to make product substitutions in lieu of repair or replacement.

For warranty service contact our customer service department at 1-800-536-0366 or at [info@beurer.com](mailto:info@beurer.com) to provide a description of the problem. If the problem is deemed to be within the scope of the limited lifetime warranty, you will be asked to mail the product at your costs in its original package with proof of purchase, your name, address and phone number. If the problem is not deemed to be within the scope of the limited lifetime warranty, we will provide a quotation for repair respectively replacement and return shipping fee.

This warranty does not cover damage caused by misuse or abuse; accident; the attachment of unauthorized accessory; alteration to the product; improper installation; misapplication; lack of reasonable care with respect to the product; unauthorized repairs or modifications; improper use of power supply; old worn battery; normal wear; loss of power; dropped product; malfunction or damage of an operating part as a result of failure to comply with instructions for use or to provide manufacturer's recommended maintenance; transit damage; theft; neglect; vandalism; or environmental conditions; loss of use during the period the product is at a repair facility or otherwise awaiting parts or repair; or any other conditions whatsoever that are beyond the control of

Beurer. This warranty is void if the product is ever used in a commercial or business environment. The maximum liability of Beurer under this warranty is limited to the purchase price actually paid by the customer for the product covered by the warranty, as confirmed by proof of purchase, regardless of the amount of any other direct or indirect damage suffered by the customer.

This warranty is effective only if the product is purchased and operated in the country in which the product is purchased. A product that requires modifications or adaptation to enable it to operate in any other country than the country for which it was designed, manufactured, approved and/or authorized, or repair of products damaged by these modifications is not covered under this warranty.

**THE WARRANTY PROVIDED HEREIN SHALL BE THE SOLE AND EXCLUSIVE WARRANTY. ANY IMPLIED WARRANTIES, OBLIGATIONS, OR LIABILITIES, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE LIMITED IN DURATION TO THE DURATION OF THIS APPLICABLE WRITTEN WARRANTY.** Some states do not allow limitations on how long an implied warranty lasts, so the above limitations may not apply to you.

**IN NO EVENT SHALL BEURER BE LIABLE FOR ANY SPECIAL, INCIDENTAL, INDIRECT OR CONSEQUENTIAL DAMAGES FOR BREACH OF THIS OR ANY OTHER WARRANTY, EXPRESS, IMPLIED OR ANY OTHER THEORY OF LIABILITY, WHATSOEVER.** Some states do not allow the exclusion or limitation of special, incidental, or consequential damages, so the above limitation may not apply to you.

Beurer does not authorize anyone, including, but not limited to, retailers, the subsequent consumer purchaser of the product from a retailer or remote purchaser, to obligate Beurer in any way beyond the terms set forth herein.

This warranty does not extend to the purchase of opened, used, repaired, repackaged and/or resealed products, including but not limited to sale of such products on Internet auction sites and/or products by surplus or bulk resellers. Any and all warranties or guarantees shall immediately cease and terminate in connection with any products or parts thereof which are repaired, replaced, altered, or modified, without the prior explicitly written consent of Beurer.

This warranty gives you specific legal rights, and you may also have other rights which may vary from state to state.

**For more information regarding our product line in the USA,  
please visit: [www.beurer.com](http://www.beurer.com)**

Beurer North America LP, Hallandale Beach, FL 33009, [www.beurer.com](http://www.beurer.com)  
Questions or comments? **Call toll-free 1-800-536-0366 or contact [info@beurer.com](mailto:info@beurer.com)**  
Made in China.

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